



# Hepacide Quat<sup>®</sup> II

**virucidal disinfectant**

Hepacide quat II is a ready-to-use, quaternary-based virucidal disinfectant spray. It is non-viscous and clear in color, and fragrance-free. Users will appreciate the following product features, advantages and benefits: effective against hepatitis, saves time and labor, user friendly and surface safe.



Providing Quality Maintenance and Chemical Specialty Solutions with Unparalleled Customer Service



Spartan Chemical Company, Inc.  
 1110 Spartan Drive  
 Maumee, OH 43537  
 1-800-537-8990  
 www.spartanchemical.com

**EFFECTIVE IN THE MEDICAL FIELD:**

HCV (Hepatitis C Virus) and HBV (Hepatitis B virus) are serious concerns in the health care environment due to its highly infectious nature and the fact that chronic infection can lead to cirrhosis and cancer of the liver. Hard water performance is not a concern since Hepacide Quat II is ready-to-use and formulated with soft water.

**NO WASTE/SAVES TIME AND LABOR:**

Hepacide Quat II does not require mixing. Since it has been registered by the EPA in its ready-to-use form, it is not necessary to make up fresh solutions each week. Unlike some concentrated disinfectants that lose their potency once mixed and must be replaced on a regular, sometimes weekly basis, Hepacide Quat II is a stable formula for up to 12 months.

**KILLS GERMS:**

Effective against the following Bacteria:

- Staphylococcus aureus
- Pseudomonas aeruginosa
- Salmonella enterica

Effective against the following Antibiotic-Resistant Bacteria:

- Vancomycin resistant Enterococcus faecalis (VRE)
- Methicillin resistant Staphylococcus aureus (MRSA)

Effective against the following Viruses:

- HIV-1 (AIDS Virus)
- Hepatitis C Virus (HCV)
- Hepatitis B Virus (HBV)
- Herpes simplex Type 2
- Influenza A2/Hong Kong viruses
- Pandemic 2009 H1N1 influenza A virus

**AREAS FOR USE:**

- |                        |                       |
|------------------------|-----------------------|
| Hospitals              | Day Care Centers      |
| Laboratories           | Public Transportation |
| Nursing Homes          | Ambulance Services    |
| Urgent Care Facilities | Law Enforcement       |
| Doctor Offices         | (Cars and Facilities) |
| Dentist Offices        | Airplanes             |
| Veterinary Clinics     | Highway Service       |
| Dialysis Centers       | Areas                 |
| Blood Banks            | Gas Stations          |
| Hotels/Motels          | Restaurants           |
| Locker Rooms           | Delicatessens         |
| Health Spas            | Cafeterias            |
| Gyms                   | Fast Food Chains      |
| Schools                | Industry              |
|                        | All Public Places     |

**APPLICATIONS:**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**Disinfection and Cleaning:** Remove gross filth before disinfection. Apply Hepacide Quat II with a cloth, mop, sponge or coarse sprayer. Treated surfaces must remain wet for 10 minutes. Wipe off with a damp cloth or rinse with water. This product will not leave grit or soap scum.

**\*\*KILLS HCV, HBV AND HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/ objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Hepatitis C Virus (HCV), Hepatitis B Virus (HBV) and human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

**Special Instructions for Cleaning and Decontamination Against HCV, HBV and HIV-1 on Surfaces or Objects Soiled with Blood or Body Fluids:**

**Personal Protection:** Exposure to blood or body fluids should be avoided by use of barrier protection such as latex or other impervious gloves. Face shields and impervious gowns should also be worn if splashing of blood or body fluids may occur.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

**Contact Time:** Leave surfaces wet for 10 minutes.

**Disposal of Infectious Material:** Contaminated cleaning materials, blood and other body fluids should be autoclaved and/or disposed of according to local regulations for infectious waste disposal.

**Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.**

**SPECIFICATIONS:**

- pH (Concentrate)– 6.0-7.0
- Color – Clear
- Scent – Fragrance Free
- Stability:
  - a. Shelf @ 24°C/75°F – One year minimum
  - b. Freeze/Thaw – Three cycles
- Phosphate Free
- EPA Reg. No. 5741-18.



Distributed by:

PACKAGING	
125165	330-gallon tote
125155	55-gallon drum
125130	30-gallon drum
125115	15-gallon drum
125105	5-gallon pail
325103	1-quart (12 per case)



**GUARANTEE:** Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 2/13 L1251