

## Glossary

### ARMREST TYPES

#### Full Length

Provide more support for easier ingress and egress.

#### Desk Length

Allow for easy fit under most desks or tables.

### RIGGING STYLES

#### Swing-Away Footrests

Provide foot support while seated and rotate out of the way for easier transfers.

#### Elevating Legrests

Offer leg support with comfortable calf pads and allow angle-of-elevation adjustment.

### CHASSIS FINISH

#### Chrome

Resistant to corrosion, easy to clean, and sturdy.

#### Powder-Coated

Smooth finish with an extremely durable coating.

## Measuring Guidelines

A wheelchair user's health and wellness is dependent upon achieving the proper fit in their new wheelchair. An improper fit can cause posture issues, breathing problems, pressure ulcers, and discomfort. Make sure the user is seated in an upright position before taking any measurements. Once you have taken the resident's measurements, use the dimensions to determine the wheelchair that will best fit the resident's needs.

#### Seat Depth Measurement

Measure from the back of the hip to the back of the knee. Subtract one inch from this measurement to determine the best wheelchair seat depth dimension for your resident.

#### Seat Width Measurement

Measure between the widest point on the hips or thighs. Use this measurement to determine the appropriate wheelchair seat width for your resident.

#### Footrest/Legrest Length Measurement


Measure from the back of the knee to the heel of the foot. Elevating legrests are always suggested when residents are required to raise their legs regularly.

#### Seat to Floor Height Measurement

Measure the distance between the back of the knee and the heel of the foot. This measurement will determine if a dual-axle wheelchair with varying seat-to-floor heights would be appropriate for the resident. If the resident will be self-propelling, a hemi seat height is recommended.

#### Back Height Measurement

Measure from the resident's collarbone to the seat. Use this measurement to determine the appropriate wheelchair back height. Taller back heights may be required for residents who need upper trunk support.










ReliaCare

Exclusively Available at  
**HDSUPPLY**

Our exclusive ReliaCare wheelchairs provide easy, tool-free adjustability, a three-year manufacturer's warranty, and a number of value-added features such as anti-tippers, an essential component for maximizing resident safety.

See our full-line Catalog V64, pages 1634-1637, and our Healthcare Catalog V9, pages 18-21, for our complete selection of wheelchairs.



TYPE	BRAND	PRODUCT NAME	WEIGHT CAPACITY	SEAT WIDTH	SEAT DEPTH	BACK HEIGHT	ARMREST STYLES	RIGGING STYLE	CHASSIS FINISH	ADJUSTABLE HEIGHT	WARRANTY	FREE DELIVERY	PART #
MANUAL	 Exclusively Available at 	Manual Wheelchair	250 Lbs	18"	16"	18"	Desk Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680615
		Manual Wheelchair	250 Lbs	18"	16"	18"	Full Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680616
		Manual Wheelchair	250 Lbs	20"	16"	18"	Desk Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680617
		Manual Wheelchair	250 Lbs	20"	16"	18"	Full Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680618
	 Yes, you can.®	Tracer EX2	250 Lbs	18"	16"	18"	Full Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554039
		Tracer EX2	250 Lbs	18"	16"	18"	Desk Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554040
		Tracer EX2	250 Lbs	20"	16"	18"	Full Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554041
		Tracer EX2	250 Lbs	20"	16"	18"	Desk Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554042
		Tracer IV	350 Lbs	20"	18"	18"	Full Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554043
		Tracer IV	350 Lbs	20"	18"	18"	Desk Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554044
		Tracer IV	350 Lbs	22"	18"	18"	Full Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554045
		Tracer IV	350 Lbs	22"	18"	18"	Desk Length	Elevating Legrests Or Swing-Away Footrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554046
 Yes, you can.®	Tracer SX5	250 Lbs	18"	16"	33"	Desk Length	Swing-Away Footrests Or Elevating Legrests	Powder-Coated Silver Vein	Yes	5-Year Frame And Cross Braces	No, Quick Ship 3-Day	554047	
 Exclusively Available at 	Manual Wheelchair	450 Lbs	22"	18"	18"	Full Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680619	
	Manual Wheelchair	450 Lbs	24"	18"	18"	Full Length	Elevating Legrests	Powder-Coated Silver Vein	Yes	3-Year Limited	Yes, Next Day	680620	
	 Bariatric Rehab Leadership	Regency XL	600 Lbs	24, 26 And 28"	18"	19-1/2"	Full Length	Swing-Away Footrests	Powder-Coated Black	Yes	Limited Lifetime Warranty On Frame And Cross Braces	No	676420, 676421, 676422
		Regency XL 2000	750 Lbs	28, 30 And 32"	20"	19-1/2"	Full Length	Swing-Away Footrests	Powder-Coated Black	Yes	Limited Lifetime Warranty On Frame And Cross Braces	No	676423, 676424, 676425